





Nebulizer Kit & Adult Mask & Air Tube Set
NEB-NSET4-83AP
Instrucion Manual

Nebulizer Kit



Components

Adult Mask



Air Tube

Exclusive parts for Omron Compressor Nebulizer, NE-C803
These parts comply with EN ISO 27427:2019.

Intended Use

This set should be used only for NE-C803, Omron compressor nebulizer.
Follow this instruction manual and applicable OMRON device manual thoroughly before using the unit.
Keep for future reference.
The nebulizer kit and its attachments are single patient multiple use.

Manufacturer 	OMRON HEALTHCARE Co., Ltd. 53, Kunotsubo, Terado-cho, Muko, KYOTO, 617-0002 JAPAN
EU-representative <div><div>EC</div><div>REP</div></div>	OMRON HEALTHCARE EUROPE B.V. Scorpius 33, 2132 LR Hoofddorp, THE NETHERLANDS www.omron-healthcare.com
Importer in EU	
Production facility	OMRON HEALTHCARE MANUFACTURING VIETNAM CO., LTD. Binh Duong Province, VIETNAM
Asia Pacific HQ	OMRON HEALTHCARE SINGAPORE PTE LTD. www.omronhealthcare-ap.com

8500180-3D
Issue Date: 2022-03-15

Made in Vietnam

Important Safety Instructions



Warning

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

(Usage)

- If you experience any allergic reactions or other difficulties during use, stop using the device immediately and consult your doctor.
- Keep the device out of reach of unsupervised infants and children. The device may contain small pieces that can be swallowed.
- Do not store the air tube while there is moisture or medication remaining inside.
- Always dispose of any remaining medication after use, and use fresh medication each time.
- Store the parts in a clean location to avoid infection.

(Cleaning and disinfecting)

Observe the rules below when cleaning or disinfecting parts. Failure to observe these rules may result in damage, inefficient nebulization or infection. For the instructions, refer to "Cleaning and disinfecting" section.

- Clean and disinfect the nebulizer kit and mask before using them when as follows:
 - the first time after purchase.
 - if the device has not been used for a long period of time
 - if more than one person uses the same device
- Be sure to wash or wipe the parts after use, and ensure that they are thoroughly disinfected and dried, and stored in a clean location.
- Do not leave the cleaning solution in the parts.



Caution

Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient, or cause damage to the equipment or the other property.

(Usage)

- In order to prevent strangulation by the power cord or air tube, or injuries caused by mishandling, please provide close supervision when this device is used by, on, or near children or invalids.
- Make sure that the parts are attached correctly.
- Do not tilt the nebulizer kit at an angle of greater than 45 degrees in all directions or shake it while in use.
- Do not use or store the device while the air tube is creased.
- Do not add more than 10 ml of medication to the medication tank.
- Do not carry or leave the nebulizer kit while the medication tank contains medication.
- Do not leave the device unattended with infants or persons who cannot express their consent.
- Do not jab the nozzle of the medication tank with a pin or any sharp object.
- Do not leave the parts in extreme hot or cold temperature, or under direct sunlight.
- To avoid the medication residue on the face, be sure to wipe the face after removing the mask.
- Do not use a damaged nebulizer kit.
- Do not use by multiple patients; Use separate accessories for each patient.

(Cleaning and disinfecting)


- Do not use a microwave oven, dish dryer or hair dryer to dry the parts.
- Do not use an autoclave, EOG gas sterilizer or low temperature plasma sterilizer.
- When disinfecting parts by boiling, make sure that the container does not boil dry. Otherwise it may also result in a fire.

General Safety Precautions:

- Inspect the parts before using them each time, and check that there are no problems.
In particular, be sure to check the following:
 - That the vaporizer head, nozzle or air tube are not damaged.
 - That the nozzle is not blocked.
- Make sure that the air tube is securely attached to the compressor and nebulizer kit, and does not come loose. Twist the air tube plug slightly when inserting it into the air tube connectors to avoid the air tube disconnecting during use.

How to use

- 1 Make sure that the power switch of your compressor is in the off (○) position.
- 2 Connect the AC adapter to the compressor and plug the power plug into a power outlet.
Note: Do not place the device in a location where it is difficult to disconnect the power cord.
- 3 Remove the inhalation top together with vaporizer head from the medication tank.
Note: If the vaporizer head came off the inhalation top and fell into the medication tank, remove it out of the task.
- 4 Add the correct amount of prescribed medication to the medication tank.
- 5 Make sure that the vaporizer head is securely attached to the inhalation top.
Note: The sizes of the left- and right-hand attached parts of the vaporizer head are not the same.
If the vaporizer head is not attached in the correct orientation, it will not nebulize correctly.
- 6 Put the inhalation top back inside the medication tank.
 - 1) Align the protrusion on the inhalation top with the indent on the medication tank, as shown.
 - 2) Rotate the inhalation top clockwise until the arrow mark on it and the alignment mark on the medication tank are aligned, and they click in place.
- 7 Attach the mask to the nebulizer kit tightly.
Note: Make sure that these attachments are fitted in such a way that the handle is positioned away from the patient. Otherwise, medication may be sprayed on patient's skin or clothes.
- 8 Attach the air tube.
Push the air tube plugs firmly into the air tube connectors.



Follow the instruction manual for NE-C803 for the following procedure.

Cleaning and disinfecting

Cleaning

Clean the parts after each use.

- 1 Discard any remaining medication in the medication cup.
- 2 Wash the parts in warm water and mild detergent.
- 3 Rinse them thoroughly with clean hot tap water and allow to air dry in a clean place.

Disinfecting

Always disinfect the parts after the last treatment of the day. If the parts are heavily stained, replace them with new ones. It is expected to withstand 365 cycles of disinfection during its one-year durable period.

Read carefully the instruction manual for NE-C803 to select an appropriate method for disinfection to each part and use one of the methods described below:


- Use a commercially available disinfectant. Follow the instructions provided by the disinfectant manufacturer.

Make sure to completely remove any remaining disinfectant.

- Disinfect the parts in boiling water.

- 1 Place the parts in a container filled with sufficient water and boil it.
- 2 Allow them to sit in the boiling water at least 15 minutes.
- 3 Remove them from the container and air dry in a clean place.

DO NOT disinfect the Air tube and Adult mask in boiling water.



Removing condensation from the air tube

If there is moisture or liquid remaining in the air tube, be sure to follow the procedure below to remove any moisture from within the air tube.

- 1 Make sure that the air tube is connected to the air connector on the compressor.
- 2 Remove the air tube from the nebulizer kit.
- 3 Turn on the compressor and pump air through the air tube to expel the moisture.

Technical data

Name	Nebulizer Kit & Adult Mask & Air Tube Set
Model	NEB-NSET4-83AP
Durable Period	1 year*
Operating conditions	+10°C to +40°C / 30% to 85%RH (non-condensing) / 700 to 1060hPa
Storage and Transport conditions	-20°C to +60°C / 10% to 95%RH (non-condensing) / 700 to 1060hPa
Contents	Nebulizer Kit (Inhalation Top, Medication Tank), Adult Mask, Air Tube, Instruction Manual

* The durable period is shown above, provided the product is used to nebulize 2 ml of medication twice a day for 7.5 minutes each time at room temperature (23°C).
It may vary depending on usage environment. Frequent use of the product may shorten the durable period.

CE 0197

Notes:

- Subject to technical modification without prior notice.
- Please report to the manufacturer and the competent authority of the Member State in which you are established about any serious incident that has occurred in relation to this device.

Description of symbols			
	Refer to the instruction manual		Temperature limitation
	Humidity limitation		Atmospheric pressure limitation
	Date of manufacture		Medical device
	LOT number		Reference catalogue number
	CE marking		